

Integration of O-18 labeling and solution isoelectric focusing in a shotgun analysis of mitochondrial proteins



Jinshan Wang

Nov. 27, 2007

Introduction

- ❁ Proteome: the time- and cell-specific protein complement of the genome
- ❁ Proteomics: large-scale study of proteins concurrently and rapidly
 - ❖ Bridge the gap between the genomic information and biological function
 - ❖ Search for disease biomarker, develop better diagnostic tests and identify therapeutic targets



Challenges

Sample preparation

- ❁ The number of proteins in a human cancer cell:
>100,000
- ❁ Dynamic range: as high as 10^{12}
- ❁ Specific classes of proteins, e.g. high pI, MW or hydrophobicity


Quantitation

- ❁ Facile methods for comparative proteomics

Proteomic Analysis

Protein preparation

Fractionation:
e.g. subcellular fractions




Protein separation and quantitation

1. Separation:

- Gel electrophoresis
- In solution separations:
 - ~ liquid chromatography
 - ~ electrophoresis
 - ~ multidimensional separation

2. Quantitation

Gel-base densitometry



Chemical labeling



Enzyme-catalyzed isotope labeling (O-18)

Protein analysis

Mass spectrometry
&
Bioinformatics

The Shotgun Approach

- ❖ The protein mixture is digested into peptides with proteases
- ❖ The peptides are separated and resolved
- ❁ Easier to be handled and separated
- ❁ Easier to fragment in the MS
- ❁ Challenging: Assemble peptides into proteins : robust separation, computational analysis and interpretation of the data

Objective

To perform a comparison study of the mitochondrial proteins between MCF-7 drug susceptible and drug resistant cell lines


- ❁ Drug resistance is the most challenging problem in chemotherapy treatment of cancer cells
- ❁ Proteomics approach for study of mechanism of drug resistance
 - ❖ increase in drug efflux (P-gp)
 - ❖ block apoptosis
- ❁ Why mitochondria?
 - ❖ central to apoptosis
 - ❖ energy factory

Outline

- ❁ Develop a reproducible method to extract soluble mitochondrial proteins from MCF-7 cells
- ❁ Combine solution isoelectric focusing (solution IEF) and reversed-phase HPLC-MS for shotgun analysis of mitochondrial peptides
- ❁ Integration of O-18 labeling with this shotgun strategy to detect changes in protein abundances between drug susceptible and mitoxantrone resistant MCF-7 cancer cells
- ❁ Evaluate this result with reverse O-18 labeling
- ❁ Consider functions of those altered proteins in drug resistance


Methods

Preparation of mitochondrial peptides



Two-dimensional separation strategy

1st dimension separation: solution IEF



The IEF device is made of multiple chambers separated by a series of microporous membrane with defined pH values.


2nd dimension separation: reversed-phase μ LC-MS/MS

NanoLC was performed using C-18 column in a linear gradient mode

Mass Spectrometry and Protein Identification

Tandem Mass Spectrometry (MS/MS)

- ❖ Mass spectrometry of fragment ions of peptides
- ❖ Very confident protein identifications







Results

Evaluation of solution IEF


? Each chamber contains different peptides




Fraction1
pH below 3




Fraction 4
pH 6.5-8




Fraction2
pH 3-5



Fraction 5
pH 8-11




Fraction3
pH 5-6.5




Fraction 6
pH above 11

? Redundant peptide distribution




? Reproducibility




UV (214nm) HPLC chromatogram of a single fraction (pH3-5) from three solution IEF separations

? Distribution of identified peptides and proteins by organelles




637 peptides were identified, of which 446 peptides originated from mitochondria




Subcellular distribution of the 278 proteins

? Distribution of pI and MW of identified mitochondrial proteins




?Comparison of identified proteins in 2D gel and this study




Quantitative Comparison

Enzymatic ^{18}O incorporation





Relative quantitation with O-18 labeling




- ❁ Global labeling strategy (exception: peptide from C-terminus of the protein)
- ❁ Compatible with analysis of proteins from tissue and other limited samples
- ❁ The by-product is water and the immobilized catalyst can be readily removed.
- ❁ Compatible with peptide separation strategies

O-18 Labeling Efficiency




Lysozyme peptide: WWCNDGR

Spectrum of $^{16}\text{O}/^{18}\text{O}$ mixture (1:1)




Quantitation of $^{18}\text{O}/^{16}\text{O}$ ratios

Ratios are calculated from mass spectral peak area





$$^{18}\text{O}/^{16}\text{O} = \{I_4 - (M_4/M_0)I_0 - M_2/M_0[I_2 - (M_2/M_0)I_0] + [I_2 - (M_2/M_0)I_0]\}/I_0$$

Integration of O-18 labeling with two-dimensional separation strategy





Unchanged abundance in drug resistant cells



Two pairs of peptides from heat shock protein 70


Altered abundance in drug resistant cells





Two pairs of peptides from elongation factor Tu

Reverse O-18 labeling



- ❖ **Forward labeling:** Peptides from drug susceptible cell line were labeled with O-18 water.
- ❖ **Reverse labeling:** Peptides from drug resistant cell line were labeled with O-18 water
- ❁ Estimate the precision and confirm the behavior of “on/off” proteins
- ❁ Help interpretation of tandem mass spectra and thus provide peptide identification




Example of quantitative comparisons based on peptide
MX: mitoxantrone-resistant cells, WT: drug-susceptible cells



Both peptides come from protein: Galectin-3 binding protein – “on/off”



Forward O-18 labeled



Reverse O-18 labeled

Galectin-3 binding protein

Proteins with altered abundances in drug resistance

Accession number	Protein name	Abundance ratio (MX/WT) in forward labeling	Abundance ratio (MX/WT) in reverse labeling
Q9HCCO	Methylcrotonoyl-CoA carboxylase beta chain	2.0±0.6	1.9±0.1
P49411	Elongation factor Tu	2.4±0.4	2.4±0.5
P35232	Prohibitin	2.9±0.3	2.4±0.1
Q92665	Mitochondrial 28S ribosomal protein S31	3.3±0.5	2.5±0.5
Q9BZZ5	Apoptosis inhibitor 5	3.3±0.1	2.4±0.2
P20674	Cytochrome C oxidase polypeptide Va	3.5±0.6	2.0±0.4
Q16698	2, 4-dienoyl-CoA reductase	3.6±0.3	2.1±0.0
Q08380	Galectin-3 binding protein precursor	Only present in MX	Only present in MX
P09622	Dihydrolipoyl dehydrogenase	0.46±0.12	0.42±0.11
Q9HAV7	GrpE protein homolog 1	0.48±0.12	0.42±0.07


Biological implications of abundance changes

Apoptosis

- ❖ **Galectin-3 binding protein (“on/off” protein, only in MX)**
 - tumor-secreted antigen in human breast cancer cells
 - galectin-3 may confer chemotherapy and apoptosis resistance
 - galectin-3 translocated to the mitochondrial membranes and inhibited cytochrome c release following a variety of apoptotic stimuli
- ❖ **Apoptosis inhibitor 5 (↑)**
 - a potent suppressor of E2F-dependent apoptosis

🌸 Fatty acid metabolism

❖ 2, 4-dienoyl-CoA reductase (↑)



- Warburg hypothesis: cancer cells consistently rely on glycolytic pathway to convert glucose to ATP (reducing ROS ?)
- fatty acid oxidation was shown to be reduced in cancer cells
- higher use of fatty acid for fuel in mitochondria in drug-resistant cells

Conclusion

- ❁ A two-dimensional solution phase separation method combined with forward and reverse O-18 labeling strategy has been established and successfully applied for protein qualitative and quantitative analysis
- ❁ The overall strategy identified ten proteins with abundances altered significantly in drug resistant cells. They are active involved in apoptosis, fatty acid metabolism, TCA cycle, etc. , and may be associated with drug resistance

Acknowledgments

- ❁ Dr. Catherine Fenselau
- ❁ Dr. Peter Gutierrez
- ❁ Dr. Nathan Edwards
- ❁ The group members
- ❁ NIH

